Alpinglow Summit Mat

User Manual

The Alpinglow Summit Mat is a versatile, portable red light therapy solution with flexible Dual LEDs, providing full-body therapy and enhancing wellness during exercise or relaxation.

Your Red Light Therapy Mat

Before starting your first session, please make sure to read this manual as it may explain important operating procedures. Notes: Each LED inside has three chips . Near Infrared 850nm is invisible to the human eyes.

This package includes

- 1 * Red light Therapy Mat
- 1 * Power Cord
- 1 * Adaptor
- 1 * Connecting Wire

The red light device uses a combo of red 660nm and Infrared 850nm LEDs which is chosen to give the ideal wavelengths found in the scientific research of the last three decades.

660nm red light

This light can be seen by the human eyes and targets skin cells.

850nm near-infrared light


This light cannot be seen by the human eyes and targets deep tissue.

Specifications


Size :	180.34 cm × 78.74 cm	Wave Length :	660nm / 850nm
Color :	Black	Brightness :	P1-P5
Power:	125W	Timer :	10-90min
Light Chips :	1280pcs	Weight:	4.75KG

Red & Near Infrared Light Therapy Mat

Super Large Red Light Therapy Mat

3 Chips in 1 Bead

How to Operate

Step 1

Connect the device to the controller

Step 2

Connect the controller with the Power adapter, and then plug in the outlet. (Don't use a power bank supply to charge the device. Cause it wildecrease the power and the intensity of the device)

Step 3

After power on, the mat should flat on the dry and smooth floor, then user need wear the anit-red light glasses, choose the mode you like and enjoy the therapy time.

Red Light Intensity P1-P5

Timer

10-90 minutes options

Power

Power On, 10Hz,40Hz,Power Off

Step 4

To choose different timer for each session, press the switch button first then press the timer button on the controller. This will add time in 10-minute increments up to 90 minutes. And press the power button again to turn off the timer function. To choose different intensities, press the switch button first then press the + button for higher intensity, or press the - button for lower intensity. And press the power button again to stop adjusting the intensities. (It is recommended to contact with your skin directly)

How to Best Use the Mat:

It is recommended to use the red light therapy mat 3–4 times per week, with each session lasting 30 minutes. After a few weeks of consistent use, you should notice positive results.

The mat offers three modes: full light normal mode, 10Hz pulse mode, and 40Hz pulse mode, as well as five energy levels (P1 to P5). We generally recommend starting with the full light normal mode and setting the energy level to P3, adjusting to P5 if more heat is desired.

For those with heat-sensitive skin, the 10Hz and 40Hz pulse modes are ideal. In these modes, the light pulses 10 or 40 times per second, creating a flashing effect that provides a gentle buffer for sensitive skin. If heat sensitivity is a concern, the 40Hz mode is particularly suitable.

How to Clean:

You must turn off the power first and make sure it is not working then clean it. DO NOT clean and immerse it in water. The red light therapy mat is not waterproof.

It only needs to be gently wiped with a dry towel. If there is fog inside the physical therapy mat and cannot see clearly, you can open the zipper and wipe it with a dry towel gently.

Benefits Backed by Science:

Research shows that the wavelengths in the spectrum of 660-850nm show optimal biological response.

- Muscle and Joints Pain
- Speed wound healing
- Enhanced blood circulation
- Weight loss
- Anti-inflammatory
- Decreases pain
- Aging support and wrinkle removal
- Increase collagen
- Hair Regrowth

Notice: Even though this technology doesn't have side effects, we suggest taking the following precautions:

- ▶ This device is not waterproof!
- ▶ Not intended for use by children and pregnant women.
- ▶ Don't stare at light directly.
- ▶ For optimal results and safety, it's best to use the mat during awake, active sessions where you can adjust settings as needed.

Thank you for investing in your health and well-being with Alpinglow Summit Mat.

www.alpinglow.ch

[alpinglow_switzerland](https://www.instagram.com/alpinglow_switzerland)

[Alpinglow](https://www.facebook.com/Alpinglow)

info@alpinglow.ch

+41 79 674 71 07

Grans Montana, Switzerland

www.alpinglow.ch

Designed and Tested in Switzerland

Disclaimer: Red light therapy panels are intended for general wellness purposes only and should not be used as a substitute for professional medical advice or treatment. Individual results may vary. Consult a healthcare professional before use.

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